## PATIENT GUIDE



# **IV Push Infusion**

### The Physician has prescribed you to administer \_\_\_\_\_

#### **Gather Supplies**

- Medication Medication(s) should always be infused at room temperature. If the medication requires refrigeration, remove it from the refrigerator 1-2 hours prior to administration.
- \_\_ Prefilled Sodium Chloride Syringes
- \_\_ Prefilled Heparin Syringe \_\_\_u/mL, if ordered
- Alcohol Pads

#### **Preparing and Infusing Medication**

- 1. After gathering the supplies above, prepare a clean, flat surface and wash hands.
- 2. Check the medication ball to ensure:
  - The patient name is correct,
  - The medication and dose are correct,
  - Expiration or Use By date has not passed,
  - Inspect syringe for leaks, cracks, particulate matter and clarity of medication/solution.
- 3. Remove air bubbles from syringes. Refer to Patient Infusion Welcome Packet for instructions.
- 4. Scrub the hub of the needleless connector on the IV line for 15 seconds with an alcohol pad. Allow to air dry.
- 5. Attach Prefilled Sodium Chloride Syringe by pushing and twisting to the right.
- 6. Unclamp the clamp on the IV line and slowly push the end of the syringe plunger using a "push-pause" method.
- 7. Remove the empty Sodium Chloride Syringe from the IV line.
- 8. Scrub the hub of the needleless connector on the line for 15 seconds with an alcohol pad. Allow to air dry.
- 9. Remove cap from medication syringe and attach medication syringe to the IV line by pushing and twisting to the right.
- 10. Slowly push \_\_\_\_\_ lines, then wait \_\_\_\_\_ minutes. Repeat process until syringe is empty.



#### **When Infusion Completed**

- 1. Remove empty medication syringe from the IV line.
- 2. Scrub the hub of the needleless connector on the line for 15 seconds with an alcohol pad. Allow to air dry.
- 3. Attach Prefilled Sodium Chloride Syringe by pushing and twisting to the right and slowly push the end of the syringe plunger using a "push-pause" method.
- 4. Remove the empty Sodium Chloride Syringe from the IV line.
- 5. Scrub the hub of the needleless connector on the line for 15 seconds with an alcohol pad. Allow to air dry.
- 6. If ordered, attach Prefilled Heparin Syringe by pushing and twisting to the right and slowly push the end of the syringe plunger.
- 7. Re-clamp clamp on the IV line.
- 8. Remove the empty Heparin Syringe from the IV line.
- 9. Dispose of all supplies including syringes in regular garbage.

#### **Questions, Concerns or Problems**

Should there be any questions, concerns or problems related to medication administration or the IV line contact the home care agency \_\_\_\_\_\_ or Upstate HomeCare.

